MercuryMedical **EZflow®**
Continuous Nebulizer (en)
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercury-med.com

INTENDED USE:
To be used under proper medical supervision to deliver medicinal aerosol and diagnostic formulations.

- CAUTION:**
- Federal law (USA) restricts this device to sale by or on the order of a physician.
 - Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
 - The medicinal product which the device is intended to administer must be approved and used under the specific direction of a physician.
 - While using supplemental oxygen, do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable chemicals.
 - Performance information may not apply to drugs supplied in suspension or high viscosity form.
 - Disposal of this device should be handled in accordance with hospital protocol.
 - This device is assembled with PVC components which include low levels of a phthalate as indicated. Consideration should be used when treating children or pregnant or nursing women.
 - Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


INSTRUCTIONS FOR CONTINUOUS NEBULIZATION:

1. Remove EZflow® from packaging.
2. Connect EZflow® to a flowmeter using oxygen supply tubing.

- CAUTION:**
- Be sure all connections are secure. Check for large tubing leaks by turning flowmeter to maximum and verifying that the flowmeter does not indicate more than 6 L/min. Check the orifice of the nebulizer and clear any obstructions.
 - 3. Pour prescribed medication into nebulizer reservoir.
 - 4. Adjust flowmeter to 3 L/min. Nominal output is 6 mL/hour.

- CAUTION:**
- The listed output is nominal value only, actual output may vary depending on device, accuracy of flowmeter, and the existence of tubing connection leaks. Follow up output checks are recommended and flowmeter adjustments may be necessary.
 - 5. Connect aerosol mask directly to the outlet of the EZflow®. Place aerosol mask on patient and ensure a good mask seal. If additional oxygen is desired, use nasal cannula under the aerosol mask.
 - 6. Check liquid consumption rate (output) every hour.
 - 7. Increase or decrease the nebulizer flow rate by small increments to obtain desired actual output.

EZFLOW® CONTINUOUS NEBULIZER SETUP:

SPECIFICATIONS:

Sample Guidelines for Preparing Rx Dosing:

5mg (1cc/Hr.)				10mg (2cc/Hr.)			
Time	Output @ 3Lpm	Albuterol in cc's	Saline Mixture	Time	Output @ 3Lpm	Albuterol in cc's	Saline Mixture
1 Hr	6cc	1cc	5cc	1 Hr	6cc	2cc	4cc
2 Hr	12cc	2cc	10cc	2 Hr	12cc	4cc	8cc
3 Hr	18cc	3cc	15cc	3 Hr	18cc	6cc	12cc
4 Hr	24cc	4cc	20cc	4 Hr	24cc	8cc	16cc

15mg (3cc/Hr.)				20mg (4cc/Hr.)			
Time	Output @ 3Lpm	Albuterol in cc's	Saline Mixture	Time	Output @ 3Lpm	Albuterol in cc's	Saline Mixture
1 Hr	6cc	3cc	3cc	1 Hr	6cc	4cc	2cc
2 Hr	12cc	6cc	6cc	2 Hr	12cc	8cc	4cc
3 Hr	18cc	9cc	9cc	3 Hr	18cc	12cc	6cc
4 Hr	24cc	12cc	12cc	4 Hr	24cc	16cc	8cc

EXAMPLE:
To deliver 10 mg/hour of 5 mg/mL medication;
4 Hour Duration: Place 8 mL of medication + 16 mL of Saline into nebulizer reservoir and run for 4 hours.

NOTE:
Unit dose solution (2.5 mg per 3 mL) delivers 5 mg/hour
Output of 6 mL/hour @ nominal flow of 3 L/min

MercuryMedical **EZflow®**
Kontinuerlig forstøver (da)
11300 - 49th Street North
Clearwater, Florida 33762-4807 USA
www.mercury-med.com

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