
**INDICATIONS FOR USE**  
The Mercury Medical Hyperinflation Bag System is a pulmonary-assist device intended to provide controlled or assisted ventilation to patients.

- CAUTION**  
Federal (USA) law restricts this device to sale by or on the order of a physician.
- This product may only be used by medical personnel well trained in techniques of positive pressure ventilation with a hyperinflation bag system.
  - New users must be trained in appropriate ventilation if a hyperinflation system is not immediately available or cannot be used effectively. (See your department manual for accepted positive pressure ventilation procedures.)
  - While using supplemental oxygen, do not allow smoking or use unit near sparking equipment, open flame, oil or other flammable materials.
  - Should not be used in toxic or hazardous atmospheres.
  - Verify proper function of the hyperinflation bag system and verify clear patient airway by monitoring for the following:
    - That the patient is being ventilated, as indicated by rise and fall of chest.
    - That the blue knurled valve rotates and the mask and patient valve are free from obstruction.
    - That clear valve obstructions, shake briskly.
    - Product should be constantly monitored while in use.
  - Gas source must be attached and flowing for hyperinflation bag system to operate properly.
    - That the breathing bag is securely and properly attached to the valve body.
  - This product should be used in conjunction with a manometer pressure gauge to determine appropriate delivery pressure. Pressure may be monitored with Mercury Medical manometer Pressure Gauge that may be incorporated on the patient valve. Units that do not feature the built-in manometer MUST be monitored with a line attached from the pressure port to a manometer pressure gauge.
  - When using the hyperinflation system, it is recommended that clinicians do not exceed 10 liters of flow per minute or 40 cm H<sub>2</sub>O pressure unless the user indicates that higher pressures are required to effectively ventilate the patient based on clinical experience and judgement.
  - Some hyperinflation units are equipped with a pressure-limiting device which opens at the pressure of 40 cm H<sub>2</sub>O. These units are equipped with a lock clip. Placing the lock clip on the pressure relief valve will prevent pressure limiting valve. Also, an abrupt, high-volume inspiratory delivery may cause the unit to exceed this level.
  - Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
  - Disposal of this device should be handled in accordance with hospital protocol.
  - This device is assembled with PVC components which include low levels of a phthalate as indicated. Consideration should be used when treating children or pregnant or nursing women.

**Figure #1:** Hyperinflation Bag with Pressure Manometer, Front View  
**Figure #2:** Hyperinflation Bag with Pressure Manometer and Optional Pressure Relief Valve, Side View

**PRE-USE CHECK (Test the Device before Use):**

- Connect oxygen tubing to flowmeter.
- Adjust the flow on the flowmeter between 5 and 10 LPM.
- Connect an infant/neonatal test lung to the patient port (If test lung is not available occlude patient port using gloved hand or thumb).
- Check the inspiratory pressure, by squeezing the bag targeting approximately 20 cm H<sub>2</sub>O on manometer. Release the breathing bag to allow it to refill (Units that do not feature the built-in manometer MUST be monitored with a line attached from pressure port to manometer pressure gauge).
- Adjust variable blue PEEP knob to desired value as indicated on the manometer. Adjusting input flow will affect PEEP. Device can achieve high PEEP, always verify with the manometer.
- Some Hyperinflation bags may come with a pressure-limiting (pop-off) device to limit pressure to approximately 40 cm H<sub>2</sub>O, thereby reducing the risk of stomach distension. If higher inspiration pressures might be necessary, test the pop-off by engaging the lock-clip. Squeeze bag as mentioned in Step #4 but target inspiratory pressure over 40 cm H<sub>2</sub>O.

**DIRECTIONS FOR USE**

- If a ventilation bag is not included, attach a ventilation bag to the patient valve.
- Oxygen supply tubing is attached to the oxygen inlet port. Attach opposite end of oxygen tubing to a flow regulated oxygen source.
- Units without the optional manometer, attach pressure tubing to the pressure port and to a manometer pressure gauge. If the unit has a built-in single-patient-use manometer, pressure tubing is not necessary.
- Set the blue knurled valve to the full open position (counterclockwise).
- Adjust the flowmeter to the desired level while ensuring that oxygen flow is adequate for ventilatory requirements.
- Attach the patient connection to the appropriate device for delivering gas to the patient.
- Squeeze the breathing bag to deliver a breath. Release the breathing bag to allow it to refill. Observe that the chest rise and fall in accordance to the desired level while ensuring that oxygen flow is adequate for ventilatory requirements.
- Adjust the patient connection to the appropriate device for delivering gas to the patient.
- Squeeze the breathing bag to deliver a breath. Release the breathing bag to allow it to refill. Observe that the chest rise and fall in accordance to the desired level while ensuring that oxygen flow is adequate for ventilatory requirements.
- Some Hyperinflation units may feature a patient valve with a special pressure-limiting device mounted on the upper valve housing. If inspiration results with pulmonary resistance, venting will occur, limiting pressure to 40cm +/-10, thereby reducing the risk of stomach distension. A hissing sound can be heard as the device opens.

NOTE: When higher inspiration pressures are necessary, the pressure limiting device can be closed with the tip of the index finger or by engaging the lock clip while squeezing the bag.

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