
(en)
LiteSaver[®]
Manometer
(Time Sequence - 10 Breaths Per Minute)

INTENDED USE:
The intended use of the LiteSaver Manometer is to provide a method to monitor airway pressure and provide a time sequence for adult patients in respiratory or cardiac arrest only with an advanced airway. (Based on 2020 AHA - American Heart Association Guidelines).

ACCESSORIES:
Mercury Medical In-Line Pressure Manometer Adapter (left or right)
Mercury Medical Pressure Manometer Stand and Tubing


WARNINGS:
1. Do not connect to pressure sources exceeding the limits of the manometer.
2. The use of the in-line adapters will increase airway dead space and their use is not recommended if the manometer can be connected directly to a built in manometer port.

3. When using the in-line adapter on resuscitators without built in manometer ports, the Mercury Medical Pressure Manometer must be placed between the patient valve and the patient.

4. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


DIRECTIONS FOR TIMING LIGHT:
1. Pull red tab to activate timing light.
2. Timing light will appear every six (6) seconds, prompting user to squeeze bag and deliver 10 breaths per minute.

CAUTIONS:
1. Single use only. Do not sterilize or re-use to reduce the risk of possible transmission of infection.
2. During use, the patient's condition must be continuously monitored.
3. Federal (USA) law restricts this device to sale by or on the order of a physician.
4. Time sequence for adults only.
5. Do not dispose of batteries in trash. Dispose of with other hazardous waste items.
6. If light does not flash, or does not flash at appropriate rate continue resuscitation following routine protocol.
7. Reuse of this device may pose a risk of cross-contamination and the device may not perform as intended.
8. MR unsafe items should not enter the MRI scanner room.

REGULATORY INFORMATION:
Disposal of Waste Electrical and Electronic Equipment (WEEE)
The Waste Electrical and Electronic Equipment (WEEE) Directive aims to minimize the impact of electrical and electronic goods on the environment, by increasing re-use and recycling and by reducing the amount of WEEE going to landfill. The symbol on product or its packaging signifies that this product must be disposed separately from ordinary household wastes at its end of life. Please kindly be aware that this is your responsibility to dispose electronic equipment at recycling centers so as to help conserve natural resources. Each country in the European Union should have its collection centers for electrical and electronic equipment recycling. For information about your recycling drop off area, please contact your related electrical and electronic equipment waste management authority.

SPECIFICATIONS:
Manometer Range: 0-60 cm H₂O (83-200-0103 & 83-200-0139)
0-40 cm H₂O (83-200-0315)
Manometer Variance: ±3 cm H₂O up to 15 cm H₂O
±5 cm H₂O over 15 cm H₂O
Manometer Port Size: 4.25 mm nominal
Ventilation Rate: 10 breaths per minute
L x W x H (including pull tab): 64 mm x 33 mm x 20 mm (2.51 in x 1.29 in x .78 in)
Weight: 4.3 g
Light Cadence Time: 8 hrs
Brightness of LED: 900 Millicandela (MCD)
Battery Shelf Life: 2 years
Storage Temperature: -20°C ~ 60°C (-4°F ~ 140°F) per battery performance data provided by manufacturer
Operating Temperature: 0°C ~ 60°C (32°F ~ 140°F)
Humidity Range: 45% ~ 85%

Specifications May Vary Depending on Device and Test Methods

Patent # US 8,522,618 B1; US 9,402,781 B2; US 10,098,809 B2; US D 748,773S; Other Patents Pending

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