


# ECG Patient Cables-Shielded



## **Cable Shielding**

All Bioconnect Unitized Patient Cables are constructed using fully shielded trunk cable and lead wires designed to exacting specifications with medical grade Santoprene® (TPE) jackets, copper center conductors, PE insulation and dual shielding. These low noise cables and wires minimize triboelectric noise and electrical interference. Cable yokes can be made with inductors (HF Chokes) or resistors if required. Where monitors are subject to EMI or ESU interference, cables may need to be equipped with ferrite beading.

## **Sheet Clip**

Anchors cable in place while in use.

## **Flex Reliefs**

Robust, Santoprene® flex reliefs are located at yoke, lead junctions and device ends of all unitized cables ensuring long life and reducing continuity disruption.

## **Traceability**

All Patient Cables are clearly marked with part numbers and lot control information.

## **Yoke Marking**

To facilitate lead attachment, cable yokes are labeled to clearly show channel designations and color codes.

## **Lead Wires**

Cable, yoke and lead wires are a single, molded unit which prevents accidental lead disconnection. Unitized cables are optimized for use with defibrillators, in ERs, in ICUs and Mobile Emergency Response Systems. The low profile, compact yoke is available in three and five lead models with snap and pinch lead terminations.

**Santoprene®** (TPE is the latex-free material chosen for use in the cables and leads. Yokes over-molds and flex reliefs are fabricated using the same TPE used in the wire and cable jackets to provide the most solid bond between cable components. This increases product life and patient comfort and safety. Santoprene® can be sterilized using most common methods: steam autoclave, gamma irradiation, and EtO. It is impervious to most cleaning solutions and solvents.



## **3 Lead**

Part #	Description	Connector #*	Resistors
#DB-3306	Patient Cable, GE Marquette	004	with 1 K
#DB-3316	Patient Cable, GE Marquette	004	without
#DB-3340	Patient Cable, AAMI 6 Pin	001	with 1 K
#DB-3351	Patient Cable, AAMI 6 Pin	001	without
#DB-3385	Patient Cable, Phillips (HP-Agilent), Fixed	002	with 1 K
#DB-3396	Patient Cable, Phillips (17 pin ECG), Fixed	003	with 1 K

## **5 Lead**

Part #	Description	Connector #*	Resistors
#DB-3506	Patient Cable, GE Marquette	004	with 1 K
#DB-3516	Patient Cable, GE Marquette	004	without
#DB-3540	Patient Cable, AAMI 6 Pin	001	with 1 K
#DB-3551	Patient Cable, AAMI 6 Pin	001	without
#DB-3585	Patient Cable, Phillips (HP-Agilent), Fixed	002	with 1 K
#DB-3596	Patient Cable, Phillips (17 pin ECG), Fixed	003	with 1 K

## **Device Connectors**

Bioconnect stocks patient cables that are compatible with an extensive list of medical monitoring devices. All device connectors are over-molded with medical grade Santoprene® and feature flexible strain reliefs for durability and ease of use.

\*See guide at end of cross-reference

# ECG Patient Cables-Unshielded



## **Shielding**

All Bioconnect Patient Cables are constructed using cable designed to exacting specifications with medical grade Santoprene® (TPE) jacket, copper center conductors, PE insulation and dual shielding. This low noise cables minimizes triboelectric noise and electrical interference. Cable yokes can be made with inductors (HF Chokes) or resistors if required. Where monitors are subject to EMI or ESU interference, cables may need to be equipped with ferrite beading.

## **Sheet Clip**

Anchors cable in place while in use.

## **Flex Reliefs**

Robust, Santoprene® flex reliefs are located at yoke, and device ends of all cables ensuring long life and reducing continuity disruption.

## **Traceability**

All Patient Cables are clearly marked with part numbers and lot control information.

## **Yoke Marking**

To facilitate lead attachment, cable yokes are labeled to clearly show channel designations and color codes.

**Santoprene®** (TPE is the latex-free material chosen for use in the cables and leads. Yokes over-molds and flex reliefs are fabricated using the same TPE used in the wire and cable jackets to provide the most solid bond between cable components. This increases product life and patient comfort and safety. Santoprene® can be sterilized using most common methods: steam autoclave, gamma irradiation, and EtO. It is impervious to most cleaning solutions and solvents.

## **Safety Interface**

This series of three and five lead Bioconnect Patient cables are designed for use with unshielded safety DIN 42-802 lead wires.



## **3 Lead**

Part #	Description	Connector #*	Resistors
#DB-1306	Patient Cable, GE Marquette	004	with 1 K
#DB-1316	Patient Cable, GE Marquette	004	without
#DB-1340	Patient Cable, AAMI 6 Pin	001	with 1 K
#DB-1351	Patient Cable, AAMI 6 Pin	001	without
#DB-1385	Patient Cable, Phillips (HP-Agilent)	002	with 1 K
#DB-1396	Patient Cable, Phillips (17 pin ECG)	003	with 1 K

## **5 Lead**

Part #	Description	Connector #*	Resistors
#DB-1506	Patient Cable, GE Marquette	004	with 1 K
#DB-1516	Patient Cable, GE Marquette	004	without
#DB-1540	Patient Cable, AAMI 6 Pin	001	with 1 K
#DB-1551	Patient Cable, AAMI 6 Pin	001	without
#DB-1585	Patient Cable, Phillips (HP-Agilent)	002	with 1 K
#DB-1596	Patient Cable, Phillips (17 pin ECG)	003	with 1 K

## **Device Connectors**

Bioconnect stocks patient cables that are compatible with an extensive list of medical monitoring devices. All device connectors are over-molded with medical grade Santoprene® and feature flexible strain reliefs for durability and ease of use.

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
<b>A D INSTRUMENTS</b>					
Mac Lab (ECG).....	001	1340	1540	3340	3540
AAEM .....	001	1340	1540	3340	3540
<b>ABBOT MEDICAL</b>					
AAXM .....	001	1340	1540	3340	3540
ABBV .....	001	1340	1540	3340	3540
ABBY Series .....	001	1340	1540	3340	3540
ABEM .....	001	1340	1540	3340	3540
ABTM .....	001	1340	1540	3340	3540
ARTM .....	001	1340	1540	3340	3540
DHS .....	001	1340	1540	3340	3540
EK-1 .....	001	1340	1540	3340	3540
EK-2 .....	001	1340	1540	3340	3540
EK-4 .....	001	1340	1540	3340	3540
EK40 .....	001	1340	1540	3340	3540
EK-40 .....	001	1340	1540	3340	3540
EK-5 .....	001	1340	1540	3340	3540
EK50 .....	001	1340	1540	3340	3540
EK-50 .....	001	1340	1540	3340	3540
EKP-1 .....	001	1340	1540	3340	3540
EKR50 .....	001	1340	1540	3340	3540
EO-8 .....	001	1340	1540	3340	3540
EX-50 .....	001	1351	1551	3351	3551
EX50 .....	001	1340	1540	3340	3540
Lifecare 5100 .....	001	1351	1551	3351	3551
Lifecare 5200 .....	001	1351	1551	3351	3551
<b>ACUSON</b>					
128XP/10C (Ultrasound) ....	001	1340	1540	3340	3540
128XP/3 (Ultrasound) .....	001	1340	1540	3340	3540
Aspen .....	001	1340	1540	3340	3540
Sequoia .....	001	1340	1540	3340	3540
<b>ADVANCED TECHNOLOGY LABS</b>					
400 .....	001	1340	1540	3340	3540
Ultramark 4 (PLUS) .....	001	1340	1540	3340	3540
Ultramark 7 .....	001	1340	1540	3340	3540
Ultramark 9 .....	001	1340	1540	3340	3540
<b>AGILENT TECHNOLOGIES</b>					
12 Pin Connector .....	002	1385	1585	3385	3585
43100A (12-Pin) .....	002	1385	1585	3385	3585
43110A (12-Pin) .....	002	1385	1585	3385	3585
43120A (12-Pin) .....	002	1385	1585	3385	3585
78352C .....	002	1385	1585	3385	3585
78354C .....	002	1385	1585	3385	3585
78834A (Neonatal Bedside) 002	002	1385	1585	3385	3585
78834C Neonatal .....	002	1385	1585	3385	3585
Agilent A1 .....	001	1340	1540	3340	3540
Codemaster .....	002	1385	1585	3385	3585
Compact M3923A .....	001	1340		3340	
Compact M3925A .....	001	1340		3340	
Compact M3929A .....	001	1340		3340	
M1001A .....	002	1385	1585	3385	3585
M1002A .....	002	1385	1585	3385	3585
M1106A .....	002	1385	1585	3385	3585
M1175A .....	002	1385	1585	3385	3585
M1176A .....	002	1385	1585	3385	3585
M1280A .....	002	1385	1585	3385	3585
M1722A .....	002	1385	1585	3385	3585
M1722B .....	002	1385	1585	3385	3585
M1723A (12-Pin) .....	002	1385	1585	3385	3585
M1723B (12-Pin) .....	002	1385	1585	3385	3585

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
Merlin System .....	002	1385	1585	3385	3585
OmniCare (CMS) .....	002	1385	1585	3385	3585
V24 .....	002	1385	1585	3385	3585
V26 .....	002	1385	1585	3385	3585
Viridia 24 .....	002	1385	1585	3385	3585
Viridia 24 C .....	002	1385	1585	3385	3585
Viridia 24 CT .....	002	1385	1585	3385	3585
Viridia M3 .....	002	1385	1585	3385	3585
Viridia M4 .....	002	1385	1585	3385	3585
<b>AIR SHIELDS</b>					
HRRM 71-1/2 .....	001	1340	1540	3340	3540
HRRM-71 .....	001	1340	1540	3340	3540
HRRM-72 .....	001	1340	1540	3340	3540
Porta-Fib I .....	001	1340	1540	3340	3540
Porta-Fib II .....	001	1340	1540	3340	3540
Porta-Fib III .....	001	1340	1540	3340	3540
System 5 (new) .....	001	1340	1540	3340	3540
System 6 (Neonatal) .....	001	1340	1540	3340	3540
<b>AMERICAN OPTICAL</b>					
1125 .....	001	1340		3340	
1700 .....	001	1340		3340	
26125 .....	001	1340		3340	
33167 .....	001	1340		3340	
5260 .....	001	1340		3340	
5330 .....	001	1340		3340	
Pulsar 4 .....	001	1340		3340	
<b>ANALOG DIGITAL INSTRUMENTS</b>					
3 Lead ECG Cable .....	001	1351	1551	3351	3551
<b>ATL</b>					
400 .....	001	1340	1540	3340	3540
Ultramark 4 (Plus) .....	001	1340	1540	3340	3540
Ultramark 7 .....	001	1340	1540	3340	3540
Ultramark 9 .....	001	1340	1540	3340	3540
<b>ATL ULTRASOUND</b>					
800 .....	001	1340	1540	3340	3540
800Plus .....	001	1340	1540	3340	3540
Apogee .....	001	1340	1540	3340	3540
CX .....	001	1340	1540	3340	3540
CX200/400 .....	001	1340	1540	3340	3540
HDI 3000 .....	001	1340	1540	3340	3540
HDI 9/2000 .....	001	1340	1540	3340	3540
Interspec Cardiology XL ....	001	1340	1540	3340	3540
Mark IV .....	001	1340	1540	3340	3540
Mark VII .....	001	1340	1540	3340	3540
Ultramark 4 .....	001	1340	1540	3340	3540
Ultramark 4 (Plus) .....	001	1340	1540	3340	3540
Ultramark 5 .....	001	1340	1540	3340	3540
Ultramark 7 .....	001	1340	1540	3340	3540
Ultramark 9 .....	001	1340	1540	3340	3540
Ultramark IX .....	001	1340	1540	3340	3540
<b>BCI</b>					
3101 .....	001	1340	1540	3340	3540
6100 .....	001	1340	1540	3340	3540
6200 .....	001	1340	1540	3340	3540
9100 .....	001	1340	1540	3340	3540
2005 .....	001	1340	1540	3340	3540
Advisor .....	001	1340	1540	3340	3540
Life Span 100 .....	001	1340	1540	3340	3540

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
<b>BECTON-DICKINSON</b>					
Guardian .....	001	1340	1540	3340	3540
Guardian 1001 .....	001	1340	1540	3340	3540
MD3A .....	001		1540		3540
N1001 .....	001	1340	1540	3340	3540
<b>BIOCHEM INTERNATIONAL</b>					
3101 .....	001	1340	1540	3340	3540
6100 .....	001	1340	1540	3340	3540
6200 .....	001	1340	1540	3340	3540
9100 .....	001	1340	1540	3340	3540
(2005) .....	001	1340	1540	3340	3540
Advisor .....	001	1340	1540	3340	3540
Life Span 100 .....	001	1340	1540	3340	3540
<b>BURDICK</b>					
208 (w/ lead select) .....	001	1340	1540	3340	3540
400 Series (5 lead) .....	001		1540		3540
4014 .....	001	1340	1540	3340	3540
4045 Defib .....	001	1340	1540	3340	3540
408 (w/ lead select) .....	001	1340	1540	3340	3540
412 (w/ lead select) .....	001	1340	1540	3340	3540
413 (w/ lead select) .....	001	1340	1540	3340	3540
414 (w/ lead select) .....	001	1340	1540	3340	3540
500 Series .....	001	1340	1540	3340	3540
511 .....	001	1340	1540	3340	3540
512 .....	001	1340	1540	3340	3540
514 .....	001	1340	1540	3340	3540
521 .....	001	1340	1540	3340	3540
522 .....	001	1340	1540	3340	3540
524 .....	001	1340	1540	3340	3540
526 .....	001	1340	1540	3340	3540
600 Series .....	001	1340	1540	3340	3540
701S .....	001	1340		3340	
702S .....	001	1340		3340	
703S .....	001	1340		3340	
900 Series - With Resistor ..	001		1540		3540
900 Series - Without Resistor	001		1551		3551
900 Series (w/ respiration) ..	001	1351	1551	3351	3551
900 Series (w/o respiration)	001	1340	1540	3340	3540
90300 .....	001		1540		3540
90303B .....	001		1540		3540
90305 .....	001		1540		3540
90308 .....	001	1340	1540	3340	3540
90308-11 .....	001	1340	1540	3340	3540
90367 Ultraview .....	003	1396	1596	3396	3596
90369 .....	003	1396	1596	3396	3596
90408 .....	001		1540		3540
90418 .....	001		1540		3540
90425 .....	001		1540		3540
90428 .....	001		1540		3540
90431 .....	001	1340	1540	3340	3540
90432 .....	001		1540		3540
90470 .....	001		1540		3540
90496 .....	003	1396	1596	3396	3596
90600 Series .....	001	1340	1540	3340	3540
90600A Series .....	001	1340	1540	3340	3540
90700 .....	001	1340	1540	3340	3540
90700 Series .....	001	1340		3340	
90701 .....	001	1340	1540	3340	3540
90721 .....	001	1340	1540	3340	3540
90724 .....	001	1340	1540	3340	3540

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
Medic 3 .....	001	1340	1540	3340	3540
Medic 4 .....	001	1340	1540	3340	3540
Medic 4 (Defib) .....	001	1340	1540	3340	3540
Medic 5 .....	001	1340	1540	3340	3540
PC Bedside .....	001	1340		3340	
PC Bed Side Series .....	001	1340		3340	
PC Express .....	001	1340		3340	
PCMS System .....	001	1340	1540	3340	3540
PCMS System 90300 Series	001	1340		3340	
TEK-208 .....	001	1340		3340	
TEK-408 .....	001	1340		3340	
TEK-412 .....	001	1340		3340	
TEK-413 .....	001	1340		3340	
TEK-413 with Lead Select ..	001	1340	1540	3340	3540
TEK-413/414 (w/ lead select)	001	1340	1540	3340	3540
TEK-414 .....	001	1340		3340	
TEK-414 with Lead Select ..	001	1340	1540	3340	3540
TEK-511 .....	001	1340	1540	3340	3540
TEK-512 .....	001	1340	1540	3340	3540
TEK-514 .....	001	1340	1540	3340	3540
TEK-521 .....	001	1340	1540	3340	3540

### CIRCAMED

ECG .....	001	1351	1551	3351	3551
-----------	-----	------	------	------	------

### CMI

BP308 .....	001	1340	1540	3340	3540
BP308E .....	001	1340		3340	
BP408 .....	001	1340	1540	3340	3540
CBM3000 .....	001	1340	1540	3340	3540
CMB-3000 .....	001	1340	1540	3340	3540

### COLIN MEDICAL INSTRUMENTS

BP308 .....	001	1340	1540	3340	3540
BP308E .....	001	1340		3340	
BP408 .....	001	1340	1540	3340	3540
CBM3000 .....	001	1340	1540	3340	3540
CMB-3000 .....	001	1340	1540	3340	3540

### COROMETRICS

1200 Defib Series .....	001	1340	1540	3340	3540
1500 Defib Series .....	001	1340	1540	3340	3540
8720 .....	001	1340	1540	3340	3540
900 Defib Series .....	001	1340	1540	3340	3540
Cinamap 8700 .....	001	1340	1540	3340	3540
Cinamap 8710 .....	001	1340	1540	3340	3540
Cinamap 8720 .....	001	1340	1540	3340	3540
Cinamap 9700 .....	001	1340	1540	3340	3540
Cinamap 9710 .....	001	1340	1540	3340	3540
Corometrics 556 .....	001	1340	1540	3340	3540
Dash 1000 .....	004	1306	1506	3306	3506
Dash 2000 .....	004	1306	1506	3306	3506
Dash 3000 .....	004	1306	1506	3306	3506
Dash 4000 .....	004	1306	1506	3306	3506
Eagle Systems .....	004	1306	1506	3306	3506
R-Wave Detector (Starcam)	001	1340	1540	3340	3540
Solar 3000 .....	004	1306	1506	3306	3506
Solar 8000 .....	004	1306	1506	3306	3506
Solar 9500 .....	004	1306	1506	3306	3506
Tram 100 .....	004	1306	1506	3306	3506
Tram 200 .....	004	1306	1506	3306	3506
Tram 300 .....	004	1306	1506	3306	3506
Tram 400 .....	004	1306	1506	3306	3506

(COROMETRICS continued on next page)

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

Bioconnect  
Connector\* 3-Lead 5-Lead  
3-Lead 5-Lead  
Shielded Shielded

Bioconnect  
Connector\* 3-Lead 5-Lead  
3-Lead 5-Lead  
Shielded Shielded

(COROMETRICS continued from previous page)

Tram 500	004	1306	1506	3306	3506
Tram 600	004	1306	1506	3306	3506
Tram 800	004	1306	1506	3306	3506
Tramscope 12	004	1306	1506	3306	3506
Tramscope 12C	004	1306	1506	3306	3506
Unity Systems	004	1306	1506	3306	3506

### CRITICARE SYSTEMS, INC.

1100	001	1340	1540	3340	3540
1508	001	1340	1540	3340	3540
504	001	1340	1540	3340	3540
504US	001	1340	1540	3340	3540
506	001	1340	1540	3340	3540
507	001	1340	1540	3340	3540
508	001	1340	1540	3340	3540
602-4 USP	001	1340	1540	3340	3540
602-A USP	001	1340	1540	3340	3540
Poet 2	001	1340	1540	3340	3540
Scholar	001	1340	1540	3340	3540

### CRITIKON INC.

1200 Defib Series	001	1340	1540	3340	3540
1500 Defib Series	001	1340	1540	3340	3540
8720	001	1340	1540	3340	3540
900 Defib Series	001	1340	1540	3340	3540
9720	001	1340	1540	3340	3540
Cinamap 8700	001	1340	1540	3340	3540
Cinamap 8710	001	1340	1540	3340	3540
Cinamap 8720	001	1340	1540	3340	3540
Cinamap 9700	001	1340	1540	3340	3540
Cinamap 9710	001	1340	1540	3340	3540
Cinamap 9720	001	1340	1540	3340	3540
Critikon 8100T	001	1340	1540	3340	3540
Corometrics 556	001	1340	1540	3340	3540
Dash 1000	004	1306	1506	3306	3506
Dash 2000	004	1306	1506	3306	3506
Dash 3000	004	1306	1506	3306	3506
Dash 4000	004	1306	1506	3306	3506
Dinamap	001	1340	1540	3340	3540
Dinamap Plus	001	1340	1540	3340	3540
Dynamap Pro	001	1340	1540	3340	3540
Eagle Systems	004	1306	1506	3306	3506
MPS	001	1340	1540	3340	3540
R-Wave Detector (Starcam)	001	1340	1540	3340	3540
Solar 3000	004	1306	1506	3306	3506
Solar 8000	004	1306	1506	3306	3506
Solar 9500	004	1306	1506	3306	3506
Tram 100	004	1306	1506	3306	3506
Tram 200	004	1306	1506	3306	3506
Tram 300	004	1306	1506	3306	3506
Tram 400	004	1306	1506	3306	3506
Tram 500	004	1306	1506	3306	3506
Tram 600	004	1306	1506	3306	3506
Tram 800	004	1306	1506	3306	3506
Tramscope 12	004	1306	1506	3306	3506
Tramscope 12C	004	1306	1506	3306	3506
Unity Systems	004	1306	1506	3306	3506

### CSI

1100	001	1340	1540	3340	3540
1508	001	1340	1540	3340	3540
504	001	1340	1540	3340	3540

504US	001	1340	1540	3340	3540
506	001	1340	1540	3340	3540
507	001	1340	1540	3340	3540
508	001	1340	1540	3340	3540
602-4 USP	001	1340	1540	3340	3540
602-A USP	001	1340	1540	3340	3540
Poet 2	001	1340	1540	3340	3540
Scholar	001	1340	1540	3340	3540

### DATAMEDIX

Guardian 1001	001	1340	1540	3340	3540
MD3A	001		1540		3540
N1001	001	1340	1540	3340	3540

### DATASCOPE

2000	001	1351	1551	3351	3551
2000 Respirator	001	1340	1540	3340	3540
2001	001	1351	1551	3351	3551
2001 Respirator	001	1340	1540	3340	3540
2002	001	1351	1551	3351	3551
2002 Respirator	001	1340	1540	3340	3540
2100	001	1351	1551	3351	3551
2100 Respirator	001	1340	1540	3340	3540
2200	001	1351	1551	3351	3551
2200 Respirator	001	1340	1540	3340	3540
3000	001	1351	1551	3351	3551
3000A	001	1351	1551	3351	3551
Ballon Pump	001	1351	1551	3351	3551
EL (Respirator)	001	1340	1540	3340	3540
Guardian	001	1340	1540	3340	3540
M/D3 w/ lead select	001	1351	1551	3351	3551
M/D3A	001	1351	1551	3351	3551
M/D4	001	1351	1551	3351	3551
Passport	001	1351	1551	3351	3551
Passport (Respirator)	001	1340	1540	3340	3540
Passport EL	001	1351	1551	3351	3551
Passport XG	001	1351	1551	3351	3551
Passport XL	001	1351	1551	3351	3551
System (SYS) 90	001	1351	1551	3351	3551
System (SYS) 90T	001	1351	1551	3351	3551
System (SYS) 97	001	1351	1551	3351	3551

### E FOR M

Guardian	001	1340	1540	3340	3540
Guardian 1001	001	1340	1540	3340	3540
MD3A	001		1540		3540
N1001	001	1340	1540	3340	3540

### ELECTRODYNE

Guardian	001	1340	1540	3340	3540
Guardian 1001	001	1340	1540	3340	3540
MD3A	001		1540		3540
N1001	001	1340	1540	3340	3540
ES-100	001	1351	1551	3351	3551

### GE MEDICAL SYSTEMS

1200 Defib Series	001	1340	1540	3340	3540
1500 Defib Series	001	1340	1540	3340	3540
556	001	1340	1540	3340	3540
8720	001	1340	1540	3340	3540
900 Defib Series	001	1340	1540	3340	3540
9720	001	1340	1540	3340	3540
Cinamap 8700	001	1340	1540	3340	3540
Cinamap 8710	001	1340	1540	3340	3540

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
Cinamap 8720	001	1340	1540	3340	3540
Cinamap 9700	001	1340	1540	3340	3540
Cinamap 9710	001	1340	1540	3340	3540
Cinamap 9720	001	1340	1540	3340	3540
Critikon 8100T	001	1340	1540	3340	3540
Corometrics 556	001	1340	1540	3340	3540
Dash 1000	004	1306	1506	3306	3506
Dash 2000	004	1306	1506	3306	3506
Dash 3000	004	1306	1506	3306	3506
Dash 4000	004	1306	1506	3306	3506
Dinamap	001	1340	1540	3340	3540
Dinamap Plus	001	1340	1540	3340	3540
Dynamap Pro	001	1340	1540	3340	3540
Eagle MAC Lab	004	1306	1506	3306	3506
Eagle Systems	004	1306	1506	3306	3506
MPS	001	1340	1540	3340	3540
R-Wave Detector (Starcam)	001	1340	1540	3340	3540
Solar 3000	004	1306	1506	3306	3506
Solar 8000	004	1306	1506	3306	3506
Solar 9500	004	1306	1506	3306	3506
Tram 100	004	1306	1506	3306	3506
Tram 200	004	1306	1506	3306	3506
Tram 300	004	1306	1506	3306	3506
Tram 400	004	1306	1506	3306	3506
Tram 500	004	1306	1506	3306	350
Tram 600	004	1306	1506	3306	3506
Tram 800	004	1306	1506	3306	3506
Tramscope 12	004	1306	1506	3306	3506
Tramscope 12C	004	1306	1506	3306	3506
Unity Systems	004	1306	1506	3306	3506

### GE TE MEDICAL

2 Channel Recorder	004	1306	1506	3306	3506
3 Channel Recorder	004	1306	1506	3306	3506

### HEWLETT PACKARD

12 Pin Connector	002	1385	1585	3385	3585
43100A (12-Pin)	002	1385	1585	3385	3585
43110A (12-Pin)	002	1385	1585	3385	3585
43120A (12-Pin)	002	1385	1585	3385	3585
78352C	002	1385	1585	3385	3585
78354C	002	1385	1585	3385	3585
78834A (Neonatal Bedside)	002	1385	1585	3385	3585
78834C Neonatal	002	1385	1585	3385	3585
Agilent A1	001	1340	1540	3340	3540
Codemaster	002	1385	1585	3385	3585
Compact M3923A	001	1340		3340	
Compact M3925A	001	1340		3340	
Compact M3929A	001	1340		3340	
M1001A	002	1385	1585	3385	3585
M1002A	002	1385	1585	3385	3585
M1106A	002	1385	1585	3385	3585
M1175A	002	1385	1585	3385	3585
M1176A	002	1385	1585	3385	3585
M1280A	002	1385	1585	3385	3585
M1722A	002	1385	1585	3385	3585
M1722B	002	1385	1585	3385	3585
M1723A (12-Pin)	002	1385	1585	3385	3585
M1723B (12-Pin)	002	1385	1585	3385	3585
Merlin System	002	1385	1585	3385	3585
OmniCare (CMS)	002	1385	1585	3385	3585
V24	002	1385	1585	3385	3585

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
V26	002	1385	1585	3385	3585
Viridia 24	002	1385	1585	3385	3585
Viridia 24 C	002	1385	1585	3385	3585
Viridia 24 CT	002	1385	1585	3385	3585
Viridia M3	002	1385	1585	3385	3585
Viridia M4	002	1385	1585	3385	3585

### HILL-ROM

Porta-Fib I	001	1340	1540	3340	3540
Porta-Fib II	001	1340	1540	3340	3540
Porta-Fib III	001	1340	1540	3340	3540
Porta-Fib III	001	1340	1540	3340	3540
System 5 (new)	001	1340	1540	3340	3540

### HONEYWELL

Guardian	001	1340	1540	3340	3540
Guardian 1001	001	1340	1540	3340	3540
MD3A	001		1540		3540
N1001	001	1340	1540	3340	3540

### ICR

208 (w/ lead select)	001	1340	1540	3340	3540
400 Series (5 lead)	001		1540		3540
4014	001	1340	1540	3340	3540
4045 Defib	001	1340	1540	3340	3540
408 (w/ lead select)	001	1340	1540	3340	3540
412 (w/ lead select)	001	1340	1540	3340	3540
413 (w/ lead select)	001	1340	1540	3340	3540
414 (w/ lead select)	001	1340	1540	3340	3540
500 Series	001	1340	1540	3340	3540
511	001	1340	1540	3340	3540
512	001	1340	1540	3340	3540
514	001	1340	1540	3340	3540
521	001	1340	1540	3340	3540
522	001	1340	1540	3340	3540
524	001	1340	1540	3340	3540
526	001	1340	1540	3340	3540
600 Series	001	1340	1540	3340	3540
701S	001	1340		3340	
702S	001	1340		3340	
703S	001	1340		3340	
900 Series - With Resistor	001		1540		3540
900 Series - Without Resistor	001		1551		3551
900 Series (w/ respiration)	001	1351	1551	3351	3551
900 Series (w/o respiration)	001	1340	1540	3340	3540
90300	001		1540		3540
90303B	001		1540		3540
90305	001		1540		3540
90308	001	1340	1540	3340	3540
90308-11	001	1340	1540	3340	3540
90367 Ultraview	003	1396	1596	3396	3596
90369	003	1396	1596	3396	3596
90408	001		1540		3540
90418	001		1540		3540
90425	001		1540		3540
90428	001		1540		3540
90431	001	1340	1540	3340	3540
90432	001		1540		3540
90470	001		1540		3540
90496	003	1396	1596	3396	3596
90600 Series	001	1340	1540	3340	3540
90600A Series	001	1340	1540	3340	3540

(ICR continued on next page)

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

Bioconnect  
Connector\* 3-Lead 5-Lead  
Shielded Shielded

Bioconnect  
Connector\* 3-Lead 5-Lead  
Shielded Shielded

(ICR continued from previous page)

90700	001	1340	1540	3340	3540
90700 Series	001	1340		3340	
90701	001	1340	1540	3340	3540
90721	001	1340	1540	3340	3540
90724	001	1340	1540	3340	3540
Medic 3	001	1340	1540	3340	3540
Medic 4	001	1340	1540	3340	3540
Medic 4 (Defib)	001	1340	1540	3340	3540
Medic 5	001	1340	1540	3340	3540
PC Bedside	001	1340		3340	
PC Bed Side Series	001	1340		3340	
PC Express	001	1340		3340	
PCMS System	001	1340	1540	3340	3540
PCMS System 90300 Series	001	1340		3340	
TEK-208	001	1340		3340	
TEK-408	001	1340		3340	
TEK-412	001	1340		3340	
TEK-413	001	1340		3340	
TEK-413 with Lead Select	001	1340	1540	3340	3540
TEK-413/414 (w/ lead select)	001	1340	1540	3340	3540
TEK-414	001	1340		3340	
TEK-414 with Lead Select	001	1340	1540	3340	3540
TEK-511	001	1340	1540	3340	3540
TEK-512	001	1340	1540	3340	3540
TEK-514	001	1340	1540	3340	3540
TEK-521	001	1340	1540	3340	3540

### INTERSPEC IREX

880	001	1351	1551	3351	3551
-----	-----	------	------	------	------

### INVIVO

400673	001	1340	1540	3340	3540
20100	001	1340	1540	3340	3540
20300	001	1340	1540	3340	3540
Escort 100	001	1340	1540	3340	3540
Escort 101	001	1340	1540	3340	3540
Escort 102	001	1340	1540	3340	3540
Escort 200	001	1340	1540	3340	3540
Escort 300	001	1340	1540	3340	3540
Escort Prism	001	1340	1540	3340	3540
Escort Prism SE	001	1340	1540	3340	3540
Prism	001	1340	1540	3340	3540

### INVIVO RESEARCH

3200	001	1340	1540	3340	3540
3202	001	1340	1540	3340	3540
3200/3202 NVS/TVS	001	1340	1540	3340	3540
M12	001	1340	1540	3340	3540
M6	001	1340	1540	3340	3540
Magnitude	001	1340	1540	3340	3540
Magnitude CV	001	1340	1540	3340	3540
Millenia	001	1351	1551	3351	3551
NVS	001	1340	1540	3340	3540
Omni Trak	001	1351	1551	3351	3551
Omni Trak 300	001	1340	1540	3340	3540
Omni Trak 3000	001	1340	1540	3340	3540
Omni Trak 3001	001	1340	1540	3340	3540
Omni Trak 3100	001	1340	1540	3340	3540
Omni Trak 3150	001	1340	1540	3340	3540
TVS	001	1340	1540	3340	3540

### IREX

880	001	1351	1551	3351	3551
-----	-----	------	------	------	------

### IVY BIOMEDICAL SYSTEMS

101	001	1351	1551	3351	3551
101 (After 1988)	001	1351	1551	3351	3551
101 (Bef 1989)	001	1351	1551	3351	3551
101 (Before 1989)	001	1340	1540	3340	3540
101 T	001	1340	1540	3340	3540
101B	001	1351	1551	3351	3551
101T	001	1351	1551	3351	3551
201	001	1351	1551	3351	3551
202	001	1351	1551	3351	3551
301	001	1351	1551	3351	3551
302	001	1351	1551	3351	3551
303	001	1351	1551	3351	3551
303A	001	1340	1540	3340	3540
303D	001	1340	1540	3340	3540
401	001	1351	1551	3351	3551
402	001	1340	1540	3340	3540
403	001	1340	1540	3340	3540
403D	001	1340	1540	3340	3540
404	001	1351	1551	3351	3551
405	001	1351	1551	3351	3551
405A	001	1340	1540	3340	3540
405D	001	1340	1540	3340	3540
405P	001	1340	1540	3340	3540
405R	001	1340	1540	3340	3540
405T	001	1340	1540	3340	3540
450	001	1340	1540	3340	3540
700	001	1340	1540	3340	3540
701	001	1340	1540	3340	3540
702	001	1340	1540	3340	3540
702A	001	1351	1551	3351	3551
703	001	1351	1551	3351	3551
703T	001	1340	1540	3340	3540
704	001	1351	1551	3351	3551
705	001	1340	1540	3340	3540
706	001	1340	1540	3340	3540
R-Guard 101 NR	001	1340	1540	3340	3540
R-Guard 101 R	001	1340	1540	3340	3540
Vital Guard 450CM	001	1340	1540	3340	3540
Vital Guard 450CN	001	1340	1540	3340	3540

### J&J

880	001	1351	1551	3351	3551
-----	-----	------	------	------	------

### LASER SYSTEM & ELECTRIC

Vista 1	001	1351	1551	3351	3551
---------	-----	------	------	------	------

### LIFE SCIENCE INSTRUMENTS

LS261	001	1340	1540	3340	3540
LS5	001	1340	1540	3340	3540

### LITTON

Guardian	001	1340	1540	3340	3540
Guardian 1001	001	1340	1540	3340	3540
MD3A	001		1540		3540
N1001	001	1340	1540	3340	3540

### LSE

Vista 1	001	1351	1551	3351	3551
LS261	001	1340	1540	3340	3540
LS5	001	1340	1540	3340	3540

### MALLINCKRODT

791018	001	1340	1540	3340	3540
4000	001	1340	1540	3340	3540

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
N-1000 .....	001	1340	1540	3340	3540
N200 .....	001	1340	1540	3340	3540
N-200 .....	001	1340	1540	3340	3540
N250 .....	001	1340	1540	3340	3540
N-250 .....	001	1340	1540	3340	3540
NPB-3900 .....	001	1340	1540	3340	3540
OEC 6120A .....	001	1340	1540	3340	3540

### MARQUETTE ELECTRONICS

1200 Defib Series .....	001	1340	1540	3340	3540
1500 Defib Series .....	001	1340	1540	3340	3540
556 .....	001	1340	1540	3340	3540
8720 .....	001	1340	1540	3340	3540
900 Defib Series .....	001	1340	1540	3340	3540
Cinamap 8700 .....	001	1340	1540	3340	3540
Cinamap 8710 .....	001	1340	1540	3340	3540
Cinamap 8720 .....	001	1340	1540	3340	3540
Cinamap 9700 .....	001	1340	1540	3340	3540
Cinamap 9710 .....	001	1340	1540	3340	3540
Corometrics 556 .....	001	1340	1540	3340	3540
Dash 1000 .....	004	1306	1506	3306	3506
Dash 2000 .....	004	1306	1506	3306	3506
Dash 3000 .....	004	1306	1506	3306	3506
Dash 4000 .....	004	1306	1506	3306	3506
Eagle MAC Lab .....	004	1306	1506	3306	3506
Eagle Systems .....	004	1306	1506	3306	3506
R-Wave Detector (Starcam) .....	001	1340	1540	3340	3540
Solar 3000 .....	004	1306	1506	3306	3506
Solar 8000 .....	004	1306	1506	3306	3506
Solar 9500 .....	004	1306	1506	3306	3506
Tram 100 .....	004	1306	1506	3306	3506
Tram 200 .....	004	1306	1506	3306	3506
Tram 300 .....	004	1306	1506	3306	3506
Tram 400 .....	004	1306	1506	3306	3506
Tram 500 .....	004	1306	1506	3306	3506
Tram 600 .....	004	1306	1506	3306	3506
Tram 800 .....	004	1306	1506	3306	3506
Tramscope 12 .....	004	1306	1506	3306	3506
Tramscope 12C .....	004	1306	1506	3306	3506
Unity Systems .....	004	1306	1506	3306	3506

### MATRIX MEDICAL

Life Defense Plus (LDP) Defib .....	001	1340	1540	3340	3540
-------------------------------------	-----	------	------	------	------

### MDE

400673 .....	001	1340	1540	3340	3540
20100 .....	001	1340	1540	3340	3540
20300 .....	001	1340	1540	3340	3540
Escort 100 .....	001	1340	1540	3340	3540
Escort 101 .....	001	1340	1540	3340	3540
Escort 102 .....	001	1340	1540	3340	3540
Escort 200 .....	001	1340	1540	3340	3540
Escort 300 .....	001	1340	1540	3340	3540
Escort Prism .....	001	1340	1540	3340	3540
Escort Prism SE .....	001	1340	1540	3340	3540
Prism .....	001	1340	1540	3340	3540

### MEDICAL DATA ELECTRONICS

400673 .....	001	1340	1540	3340	3540
20100 .....	001	1340	1540	3340	3540
20300 .....	001	1340	1540	3340	3540
Escort 100 .....	001	1340	1540	3340	3540
Escort 101 .....	001	1340	1540	3340	3540

Bioconnect  
Connector\* 3-Lead 5-Lead 3-Lead 5-Lead  
Shielded Shielded

Escort 102 .....	001	1340	1540	3340	3540
Escort 200 .....	001	1340	1540	3340	3540
Escort 300 .....	001	1340	1540	3340	3540
Escort Prism .....	001	1340	1540	3340	3540
Escort Prism SE .....	001	1340	1540	3340	3540
Prism .....	001	1340	1540	3340	3540

### MEDICAL RESEARCH LABS

360SL .....	001	1340	1540	3340	3540
450 Preamp .....	001	1340	1540	3340	3540
450SL .....	001	1340	1540	3340	3540
500PBS .....	001	1340	1540	3340	3540
521A .....	001	1340	1540	3340	3540
EMS-5 .....	001	1340	1540	3340	3540
EMS-6 .....	001	1340	1540	3340	3540
MD521 .....	001	1340	1540	3340	3540
MDR5200 .....	001	1340	1540	3340	3540
MDR521 .....	001	1340	1540	3340	3540
MDR529 .....	001	1340	1540	3340	3540
Portapak 1000 Isolated .....	001	1340	1540	3340	3540
PortaPak 80 .....	001	1340	1540	3340	3540
PortaPak 90 .....	001	1340	1540	3340	3540
Pre-Amp (Isolated) .....	001	1340	1540	3340	3540
Pre-Amp-450 .....	001	1340	1540	3340	3540

### MEDTEL

114-500 .....	001	1351	1551	3351	3551
351-126 .....	001	1351	1551	3351	3551
351-696 .....	001	1340	1540	3340	3540
HS20 .....	001	1351	1551	3351	3551
HS30 .....	001	1351	1551	3351	3551
HS40 .....	001	1351	1551	3351	3551

### MEDTRONIC PHYSIO-CONTROL

1440 .....	001	1551	3351	3551	
911ST .....	001	1551	3351	3551	
Life Pak 5 with lead select ..	001	1351	1551	3351	3551
Lifepak 10 Straight Connector ..	001	1551	3351	3551	
Lifepak 5 Straight Connector	001	1551	3351	3551	
Lifepak 6 .....	001	1351	1551	3351	3551
Lifepak 6S .....	001	1551	3351	3551	
Lifepak 7 Straight Connector	001	1551	3351	3551	
Lifepak 8 Straight Connector	001	1551	3351	3551	
Lifepak 9 Straight Connector	001	1551	3351	3551	
Lifepak 9A .....	001	1551	3351	3551	
Lifepak 9P .....	001	1551	3351	3551	
Physio 1440/R .....	001	1351	1551	3351	3551
VSM-1 .....	001	1340	1540	3340	3540
VSM-1ESF (w/ resistors) .....	001	1340	1540	3340	3540
VSM-1ESF (w/o resistors) ..	001	1351	1551	3351	3551
VSM-2 .....	001	1351	1551	3351	3551
VSM-2ESF .....	001	1351	1551	3351	3551
VSM-3 .....	001	1340	1540	3340	3540
VSM-3 (w/o resistors) .....	001	1351	1551	3351	3551
VSM-4 .....	001	1340	1540	3340	3540
VSM-5 .....	001	1340	1540	3340	3540

### MOTOROLA

Apcor Radio .....	001	1340	1540	3340	3540
-------------------	-----	------	------	------	------

### MRL

360SL .....	001	1340	1540	3340	3540
450 Preamp .....	001	1340	1540	3340	3540

(MRL continued on next page)

\*See guide at end of cross-reference



# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

Bioconnect  
Connector\* 3-Lead 5-Lead  
3-Lead 5-Lead  
Shielded Shielded

Bioconnect  
Connector\* 3-Lead 5-Lead  
3-Lead 5-Lead  
Shielded Shielded

(MRL continued from previous page)

450SL	001	1340	1540	3340	3540
500PBS	001	1340	1540	3340	3540
521A	001	1340	1540	3340	3540
EMS-5	001	1340	1540	3340	3540
EMS-6	001	1340	1540	3340	3540
MD521	001	1340	1540	3340	3540
MDR5200	001	1340	1540	3340	3540
MDR521	001	1340	1540	3340	3540
MDR529	001	1340	1540	3340	3540
Portapak 1000 Isolated	001	1340	1540	3340	3540
PortaPak 80	001	1340	1540	3340	3540
PortaPak 90	001	1340	1540	3340	3540
Pre-Amp (Isolated)	001	1340	1540	3340	3540
Pre-Amp-450	001	1340	1540	3340	3540

### NARCO

HRRM 71-1/2	001	1340	1540	3340	3540
HRRM-71	001	1340	1540	3340	3540
HRRM-72	001	1340	1540	3340	3540
Porta-Fib I	001	1340	1540	3340	3540
Porta-Fib II	001	1340	1540	3340	3540
Porta-Fib III	001	1340	1540	3340	3540
System 5 (new)	001	1340	1540	3340	3540
System 6 (Neonatal)	001	1340	1540	3340	3540

### NELCOR INC

791018	001	1340	1540	3340	3540
4000	001	1340	1540	3340	3540
N-1000	001	1340	1540	3340	3540
N200	001	1340	1540	3340	3540
N-200	001	1340	1540	3340	3540
N250	001	1340	1540	3340	3540
N-250	001	1340	1540	3340	3540
NPB-3900	001	1340	1540	3340	3540
OEC 6120A	001	1340	1540	3340	3540

### NIHON KOHDEN

BSM-8301	001	1340	1540	3340	3540
BSM-8301A	001	1340	1540	3340	3540
BSM-8302	001	1340	1540	3340	3540
BSM850 SERIES	001	1340	1540	3340	3540
BSM-8500A	001	1340	1540	3340	3540
BSM-8502	001	1340	1540	3340	3540
Cardiolife 5000	001	1340	1540	3340	3540
Cardiolife 7000	001	1340	1540	3340	3540
Life Scope 9	001	1351	1551	3351	3551
OEC-6102A (s/n>20156)	001	1340	1540	3340	3540
OEC-6105A (s/n<20226)	001	1340	1540	3340	3540
OEC-6105A (s/n>20226)	001	1340	1540	3340	3540
OEC-6120A	001	1340	1540	3340	3540
OEC-6120A	001	1340	1540	3340	3540
OEC-6120J	001	1340	1540	3340	3540
OEC-6120K	001	1340	1540	3340	3540
OEC-7102A	001	1340	1540	3340	3540
OEC-8108A	001	1340	1540	3340	3540
OEC-8120A	001	1340	1540	3340	3540
TEC-5200A (s/n>20421)	001	1340	1540	3340	3540
TEC-7100A	001	1340	1540	3340	3540
TEC-7200A	001	1340	1540	3340	3540
TEC-7300A	001	1340	1540	3340	3540

### NORTH AMERICAN DRAEGER

Cable with out Resistors	001	1351	1551	3351	3551
Cable with Resistors	001	1340	1540	3340	3540
Vitalert 1000	001	1340	1540	3340	3540
Vitalert 2000	001	1340	1540	3340	3540

### NOVAMETRIX MEDICAL SYSTEMS

902	001	1351	1551	3351	3551
903	001	1351	1551	3351	3551
910	001	1351	1551	3351	3551
Life Defense 1	001	1340	1540	3340	3540
Neonatal	001	1340	1540	3340	3540
Pace-Tech MiniMax 4000CL	001	1351		3351	
Minipack 300	001	1351		3351	
Minipack 3000	001	1351		3351	
Minipack 3100	001	1351		3351	
Vitalmax 800 Plus	001	1351		3351	
Vitalmax 2000	001	1351		3351	
Vitalmax 2100	001	1351		3351	
Vitalmax 2200	001	1351		3351	
Vitalmax 4000	001	1351		3351	
Vitalmax 4000 Modular CL	001	1351		3351	
Vitalmax 4000CL	001	1351		3351	
Vitalmax 4100	001	1351		3351	
Vitalmax 4100 CL	001	1351		3351	

### PHILIPS MEDICAL SYSTEMS

12 Pin Connector	002	1385	1585	3385	3585
43100A (12-Pin)	002	1385	1585	3385	3585
43110A (12-Pin)	002	1385	1585	3385	3585
43120A (12-Pin)	002	1385	1585	3385	3585
78352C (12-Pin)	002	1385	1585	3385	3585
78354C (12-Pin)	002	1385	1585	3385	3585
78834A					
Neonatal Bedside (12-Pin)	002	1385	1585	3385	3585
78834C Neonatal (12-Pin)	002	1385	1585	3385	3585
Agilent A1	002	1385		3385	
Codemaster	002	1385	1585	3385	3585
Compact M3923A	001	1340		3340	
Compact M3925A	001	1340		3340	
Compact M3929A	001	1340		3340	
M1001A (12-Pin)	002	1385	1585	3385	3585
M1002A (12-Pin)	002	1385	1585	3385	3585
M1106A (12-Pin)	002	1385	1585	3385	3585
M1175A (12-Pin)	002	1385	1585	3385	3585
M1176A (12-Pin)	002	1385	1585	3385	3585
M1280A (12-Pin)	002	1385	1585	3385	3585
M1722A (12-Pin)	002	1385	1585	3385	3585
M1722B (12-Pin)	002	1385	1585	3385	3585
M1723A (12-Pin)	002	1385	1585	3385	3585
M1723B (12-Pin)	002	1385	1585	3385	3585
Merlin System (12-Pin)	002	1385	1585	3385	3585
OmniCare (CMS) (12-Pin)	002	1385	1585	3385	3585
V24 (12-Pin)	002	1385	1585	3385	3585
V26 (12-Pin)	002	1385	1585	3385	3585
Viridia 24 (12-Pin)	002	1385	1585	3385	3585
Viridia 24 C (12-Pin)	002	1385	1585	3385	3585
Viridia 24 CT (12-Pin)	002	1385	1585	3385	3585
Viridia M3 (12-Pin)	002	1385	1585	3385	3585
Viridia M4 (12-Pin)	002	1385	1585	3385	3585

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

Bioconnect  
Connector\* 3-Lead 5-Lead  
3-Lead 5-Lead  
Shielded Shielded

Bioconnect  
Connector\* 3-Lead 5-Lead  
3-Lead 5-Lead  
Shielded Shielded

### PHYSIO-CONTROL

1440	001	1551	3551
911ST	001	1551	3551
Life Pak 5 with lead select	001	1351 ... 1551	3351 ... 3551
Lifepak 10 Straight Connector	001	1551	3551
Lifepak 5 Straight Connector	001	1551	3551
Lifepak 6	001	1351 ... 1551	3351 ... 3551
Lifepak 6S	001	1551	3551
Lifepak 7 Straight Connector	001	1551	3551
Lifepak 8 Straight Connector	001	1551	3551
Lifepak 9 Straight Connector	001	1551	3551
Lifepak 9A	001	1551	3551
Lifepak 9P	001	1551	3551
Physio 1440/R	001	1351 ... 1551	3351 ... 3551
VSM-1	001	1340 ... 1540	3340 ... 3540
VSM-1ESF (w/ resistors)	001	1340 ... 1540	3340 ... 3540
VSM-1ESF (w/o resistors)	001	1351 ... 1551	3351 ... 3551
VSM-2	001	1351 ... 1551	3351 ... 3551
VSM-2ESF	001	1351 ... 1551	3351 ... 3551
VSM-3	001	1340 ... 1540	3340 ... 3540
VSM-3 (w/o resistors)	001	1351 ... 1551	3351 ... 3551
VSM-4	001	1340 ... 1540	3340 ... 3540
VSM-5	001	1340 ... 1540	3340 ... 3540

### PIONEER

1125	001	1340	3340
1700	001	1340	3340
26125	001	1340	3340
33167	001	1340	3340
33167	001	1340	3340
5260	001	1340	3340
5260	001	1340	3340
5330	001	1340	3340
5330	001	1340	3340
Pulsar 4	001	1340	3340

### PPG

Guardian	001	1340 ... 1540	3340 ... 3540
Guardian 1001	001	1340 ... 1540	3340 ... 3540
MD3A	001	1540	3540
N1001	001	1340 ... 1540	3340 ... 3540

### PROTOCOL SYSTEMS

Atlas	001	1340 ... 1540	3340 ... 3540
Propack 102	001	1340 ... 1540	3340 ... 3540
Propack 103	001	1340 ... 1540	3340 ... 3540
Propack 104	001	1340 ... 1540	3340 ... 3540
Propack 106	001	1340 ... 1540	3340 ... 3540
Propaq 100 Series (Disc.)	001	1340 ... 1540	3340 ... 3540
Propaq 200 Series	001	1340 ... 1540	3340 ... 3540
Propaq Atlas			
(3,5 Ld w/ Lead Select)	001	1340 ... 1540	3340 ... 3540
Propaq CS			
(3,5 Ld w/ Lead Select)	001	1340 ... 1540	3340 ... 3540
Propaq Encore	001	1340 ... 1540	3340 ... 3540
Propaq Encore			
(3,5 Ld w/ Lead Select)	001	1340 ... 1540	3340 ... 3540
Propaq Encore 202	001	1340 ... 1540	3340 ... 3540
Propaq Encore 204	001	1340 ... 1540	3340 ... 3540
Propaq Encore 206	001	1340 ... 1540	3340 ... 3540
Ultra Smart	001	1340 ... 1540	3340 ... 3540
Ultra Smart with 90 Degrees	001	1340 ... 1540	3340 ... 3540
Vital Signs	001	1340 ... 1540	3340 ... 3540

### PURITAN BENNETT

791018	001	1340 ... 1540	3340 ... 3540
4000	001	1340 ... 1540	3340 ... 3540
N-1000	001	1340 ... 1540	3340 ... 3540
N200	001	1340 ... 1540	3340 ... 3540
N-200	001	1340 ... 1540	3340 ... 3540
N250	001	1340 ... 1540	3340 ... 3540
N-250	001	1340 ... 1540	3340 ... 3540
NPB-3900	001	1340 ... 1540	3340 ... 3540
OEC 6120A	001	1340 ... 1540	3340 ... 3540

### QUINTON INSTRUMENTS

Q-Cath	001	1340 ... 1540	3340 ... 3540
Q-Cath	001	1340 ... 1540	3340 ... 3540

### SIEMENS ELECTROMEDICAL SYSTEMS

350 Series	001	1340	3340
400-500 ECG Respirator	001	1351 ... 1551	3351 ... 3551
Sirecust 630	001	1340	3340
Sirecust DU	001	1340	3340

### SPACELABS MEDICAL

208 (w/ lead select)	001	1340 ... 1540	3340 ... 3540
400 Series (5 lead)	001	1540	3540
4014	001	1340 ... 1540	3340 ... 3540
4045 Defib	001	1340 ... 1540	3340 ... 3540
408 (w/ lead select)	001	1340 ... 1540	3340 ... 3540
412 (w/ lead select)	001	1340 ... 1540	3340 ... 3540
413 (w/ lead select)	001	1340 ... 1540	3340 ... 3540
414 (w/ lead select)	001	1340 ... 1540	3340 ... 3540
500 Series	001	1340 ... 1540	3340 ... 3540
511	001	1340 ... 1540	3340 ... 3540
512	001	1340 ... 1540	3340 ... 3540
514	001	1340 ... 1540	3340 ... 3540
521	001	1340 ... 1540	3340 ... 3540
522	001	1340 ... 1540	3340 ... 3540
524	001	1340 ... 1540	3340 ... 3540
526	001	1340 ... 1540	3340 ... 3540
600 Series	001	1340 ... 1540	3340 ... 3540
701S	001	1340	3340
702S	001	1340	3340
703S	001	1340	3340
900 Series - With Resistor	001	1540	3540
900 Series - Without Resistor	001	1551	3551
900 Series (w/ respiration)	001	1351 ... 1551	3351 ... 3551
900 Series (w/o respiration)	001	1340 ... 1540	3340 ... 3540
90300	001	1540	3540
90303B	001	1540	3540
90305	001	1540	3540
90308	001	1340 ... 1540	3340 ... 3540
90308-11	001	1340 ... 1540	3340 ... 3540
90367 Ultraview	003	1396 ... 1596	3396 ... 3596
90369	003	1396 ... 1596	3396 ... 3596
90408	001	1540	3540
90418	001	1540	3540
90425	001	1540	3540
90428	001	1540	3540
90431	001	1340 ... 1540	3340 ... 3540
90432	001	1540	3540
90470	001	1540	3540
90496	003	1396 ... 1596	3396 ... 3596
90600 Series	001	1340 ... 1540	3340 ... 3540

(SPACELABS continued on next page)

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
<i>(SPACELABS continued from previous page)</i>					
90600A Series	001	1340	1540	3340	3540
90700	001	1340	1540	3340	3540
90700 Series	001	1340		3340	
90701	001	1340	1540	3340	3540
90721	001	1340	1540	3340	3540
90724	001	1340	1540	3340	3540
Medic 3	001	1340	1540	3340	3540
Medic 4	001	1340	1540	3340	3540
Medic 4 (Defib)	001	1340	1540	3340	3540
Medic 5	001	1340	1540	3340	3540
PC Bedside	001	1340		3340	
PC Bed Side Series	001	1340		3340	
PC Express	001	1340		3340	
PCMS System	001	1340	1540	3340	3540
PCMS System 90300 Series	001	1340		3340	
TEK-208	001	1340		3340	
TEK-408	001	1340		3340	
TEK-412	001	1340		3340	
TEK-413	001	1340		3340	
TEK-413 with Lead Select	001	1340	1540	3340	3540
TEK-413/414 (w/ lead select)	001	1340	1540	3340	3540
TEK-414	001	1340		3340	
TEK-414 with Lead Select	001	1340	1540	3340	3540
TEK-511	001	1340	1540	3340	3540
TEK-512	001	1340	1540	3340	3540
TEK-514	001	1340	1540	3340	3540
TEK-521	001	1340	1540	3340	3540
<b>TEKTRONIX</b>					
208 (w/ lead select)	001	1340	1540	3340	3540
400 Series (5 lead)	001			1540	3540
4014	001	1340	1540	3340	3540
4045 Defib	001	1340	1540	3340	3540
408 (w/ lead select)	001	1340	1540	3340	3540
412 (w/ lead select)	001	1340	1540	3340	3540
413 (w/ lead select)	001	1340	1540	3340	3540
414 (w/ lead select)	001	1340	1540	3340	3540
500 Series	001	1340	1540	3340	3540
511	001	1340	1540	3340	3540
512	001	1340	1540	3340	3540
514	001	1340	1540	3340	3540
521	001	1340	1540	3340	3540
522	001	1340	1540	3340	3540
524	001	1340	1540	3340	3540
526	001	1340	1540	3340	3540
600 Series	001	1340	1540	3340	3540
701S	001	1340		3340	
702S	001	1340		3340	
703S	001	1340		3340	
900 Series - With Resistor	001		1540		3540
900 Series - Without Resistor	001		1551		3551
900 Series (w/ respiration)	001	1351	1551	3351	3551
900 Series (w/o respiration)	001	1340	1540	3340	3540
90300	001		1540		3540
90303B	001		1540		3540
90305	001		1540		3540
90308	001	1340	1540	3340	3540
90308-11	001	1340	1540	3340	3540
90367 Ultraview	003	1396	1596	3396	3596
90369	003	1396	1596	3396	3596
90408	001		1540		3540

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
90418	001		1540		3540
90425	001		1540		3540
90428	001		1540		3540
90431	001	1340	1540	3340	3540
90432	001		1540		3540
90470	001		1540		3540
90496	003	1396	1596	3396	3596
90600 Series	001	1340	1540	3340	3540
90600A Series	001	1340	1540	3340	3540
90700	001	1340	1540	3340	3540
90700 Series	001	1340		3340	
90701	001	1340	1540	3340	3540
90721	001	1340	1540	3340	3540
90724	001	1340	1540	3340	3540
Medic 3	001	1340	1540	3340	3540
Medic 4	001	1340	1540	3340	3540
Medic 4 (Defib)	001	1340	1540	3340	3540
Medic 5	001	1340	1540	3340	3540
PC Bedside	001	1340		3340	
PC Bed Side Series	001	1340		3340	
PC Express	001	1340		3340	
PCMS System	001	1340	1540	3340	3540
PCMS System 90300 Series	001	1340		3340	
TEK-208	001	1340		3340	
TEK-408	001	1340		3340	
TEK-412	001	1340		3340	
TEK-413	001	1340		3340	
TEK-413 with Lead Select	001	1340	1540	3340	3540
TEK-413/414 (w/ lead select)	001	1340	1540	3340	3540
TEK-414	001	1340		3340	
TEK-414 with Lead Select	001	1340	1540	3340	3540
TEK-511	001	1340	1540	3340	3540
TEK-512	001	1340	1540	3340	3540
TEK-514	001	1340	1540	3340	3540
TEK-521	001	1340	1540	3340	3540
<b>TEMTECH</b>					
Lifeguard 9	001	1351	1551	3351	3551
M401	001	1351	1551	3351	3551
M501	001	1351	1551	3351	3551
M503	001	1351	1551	3351	3551
<b>VIASYS</b>					
400673	001	1340	1540	3340	3540
20100	001	1340	1540	3340	3540
20300	001	1340	1540	3340	3540
Escort 100	001	1340	1540	3340	3540
Escort 101	001	1340	1540	3340	3540
Escort 102	001	1340	1540	3340	3540
Escort 200	001	1340	1540	3340	3540
Escort 300	001	1340	1540	3340	3540
Escort Prism	001	1340	1540	3340	3540
Escort Prism SE	001	1340	1540	3340	3540
Prism	001	1340	1540	3340	3540
<b>VICKERS</b>					
HRRM 71-1/2	001	1340	1540	3340	3540
HRRM-71	001	1340	1540	3340	3540
HRRM-72	001	1340	1540	3340	3540
Porta-Fib I	001	1340	1540	3340	3540
Porta-Fib II	001	1340	1540	3340	3540
Porta-Fib III	001	1340	1540	3340	3540
Porta-Fib III	001	1340	1540	3340	3540

\*See guide at end of cross-reference

# ECG Patient Cables

## Medical Device Cross-Reference

All part numbers  
begin with #DB-

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
System 5 (new) .....	001	1340	1540	3340	3540
System 6 (Neonatal) .....	001	1340	1540	3340	3540
<b>VIATEK</b>					
208 (w/ lead select) .....	001	1340	1540	3340	3540
400 Series (5 lead) .....	001		1540		3540
4014 .....	001	1340	1540	3340	3540
4045 Defib .....	001	1340	1540	3340	3540
408 (w/ lead select) .....	001	1340	1540	3340	3540
412 (w/ lead select) .....	001	1340	1540	3340	3540
413 (w/ lead select) .....	001	1340	1540	3340	3540
414 (w/ lead select) .....	001	1340	1540	3340	3540
511 .....	001	1340	1540	3340	3540
512 .....	001	1340	1540	3340	3540
514 .....	001	1340	1540	3340	3540
521 .....	001	1340	1540	3340	3540
522 .....	001	1340	1540	3340	3540
524 .....	001	1340	1540	3340	3540
526 .....	001	1340	1540	3340	3540
500 Series .....	001	1340	1540	3340	3540
600 Series .....	001	1340	1540	3340	3540
701S .....	001	1340		3340	
702S .....	001	1340		3340	
703S .....	001	1340		3340	
900 Series - With Resistor .....	001		1540		3540
900 Series - Without Resistor .....	001		1551		3551
900 Series (w/ respiration) .....	001	1351	1551	3351	3551
900 Series (w/o respiration) .....	001	1340	1540	3340	3540
90300 .....	001		1540		3540
90303B .....	001		1540		3540
90305 .....	001		1540		3540
90308 .....	001	1340	1540	3340	3540
90308-11 .....	001	1340	1540	3340	3540
90367 Ultraview .....	003	1396	1596	3396	3596
90369 .....	003	1396	1596	3396	3596
90408 .....	001		1540		3540
90418 .....	001		1540		3540
90425 .....	001		1540		3540
90428 .....	001		1540		3540
90431 .....	001	1340	1540	3340	3540
90432 .....	001		1540		3540
90470 .....	001		1540		3540
90496 .....	003	1396	1596	3396	3596
90600 Series .....	001	1340	1540	3340	3540
90600A Series .....	001	1340	1540	3340	3540
90700 .....	001	1340	1540	3340	3540
90700 Series .....	001	1340		3340	
90701 .....	001	1340	1540	3340	3540
90721 .....	001	1340	1540	3340	3540
90724 .....	001	1340	1540	3340	3540
Medic 3 .....	001	1340	1540	3340	3540
Medic 4 .....	001	1340	1540	3340	3540
Medic 4 (Defib) .....	001	1340	1540	3340	3540
Medic 5 .....	001	1340	1540	3340	3540
PC Bedside .....	001	1340		3340	
PC Bed Side Series .....	001	1340		3340	
PC Express .....	001	1340		3340	
PCMS System .....	001	1340	1540	3340	3540
PCMS System 90300 Series .....	001	1340		3340	
TEK-208 .....	001	1340		3340	
TEK-408 .....	001	1340		3340	

	Bioconnect Connector*	3-Lead	5-Lead	3-Lead Shielded	5-Lead Shielded
TEK-412 .....	001	1340		3340	
TEK-413 .....	001	1340		3340	
TEK-413 with Lead Select .....	001	1340	1540	3340	3540
TEK-413/414 (w/ lead select) .....	001	1340	1540	3340	3540
TEK-414 .....	001	1340		3340	
TEK-414 with Lead Select .....	001	1340	1540	3340	3540
TEK-511 .....	001	1340	1540	3340	3540
TEK-512 .....	001	1340	1540	3340	3540
TEK-514 .....	001	1340	1540	3340	3540
TEK-521 .....	001	1340	1540	3340	3540

### WELCH ALLYN PROTOCOL

Atlas .....	001	1340	1540	3340	3540
Propack 102 .....	001	1340	1540	3340	3540
Propack 103 .....	001	1340	1540	3340	3540
Propack 104 .....	001	1340	1540	3340	3540
Propack 106 .....	001	1340	1540	3340	3540
Propaq 100 Series (Disc.) .....	001	1340	1540	3340	3540
Propaq 200 Series .....	001	1340	1540	3340	3540
Propaq Atlas (3,5 Ld w/ Lead Select) .....	001	1340	1540	3340	3540
Propaq CS (3,5 Ld w/ Lead Select) .....	001	1340	1540	3340	3540
Propaq Encore .....	001	1340	1540	3340	3540
Propaq Encore (3,5 Ld w/ Lead Select) .....	001	1340	1540	3340	3540
Propaq Encore 202 .....	001	1340	1540	3340	3540
Propaq Encore 204 .....	001	1340	1540	3340	3540
Propaq Encore 206 .....	001	1340	1540	3340	3540
Ultra Smart .....	001	1340	1540	3340	3540
Ultra Smart with 90 Degrees .....	001	1340	1540	3340	3540
Vital Signs .....	001	1340	1540	3340	3540

### ZMI

8000-00-25 .....	001	1340	1540	3340	3540
D-90 .....	001	1340	1540	3340	3540
D900 .....	001	1340	1540	3340	3540
PD1200 .....	001	1340	1540	3340	3540
PD1400 .....	001	1340	1540	3340	3540
PD-2000 .....	001	1340	1540	3340	3540

### ZOLL MANUFACTURING (ZMI)






8000-00-25 .....	001	1340	1540	3340	3540
D-90 .....	001	1340	1540	3340	3540
D900 .....	001	1340	1540	3340	3540
PD1200 .....	001	1340	1540	3340	3540
PD1400 .....	001	1340	1540	3340	3540
PD-2000 .....	001	1340	1540	3340	3540

This cross-reference guide is based upon the best data available at the time of compilation from many different sources. Devices may change without notice. Please match your device connector to the illustrations found at the end of this cross-reference.


\*See guide at end of cross-reference

# ECG Patient Cables


## Reference Guide




1300 Series Patient Cable




1500 Series Patient Cable



DIN Style Safety Lead



3300 Series Shielded Patient Cable



3500 Series Shielded Patient Cable

These illustrations are provided to confirm that our connector is a match for your device.  
Carefully count the exact number of pins or sockets.