


Miscellaneous



DIAGNOSTIC LIGHTS with Pupil Gauge

#10-000-01 6/pk

This small, lightweight penlight clips to the pocket. It weighs less than one ounce, yet delivers a bright dependable beam. It is perfect for office exams, clinical rounds and emergencies. To activate the light, press the clip switch. The penlight will last from one to six months, depending upon usage. When the batteries are depleted, the penlight is discarded. The Mercury six-pack is packaged in a self-sealing poly bag for refrigerator storage of unused lights. Refrigeration will keep lights fresh.



FAS® Clip #69-20745

Reopenable, adjustable clip ideal for cables, hoses, etc. Adhesive backing.