



NO POSTAGE  
NECESSARY  
IF MAILED  
IN THE  
UNITED STATES



**BUSINESS REPLY MAIL**

FIRST-CLASS MAIL PERMIT NO. 412 ST. PETERSBURG, FL  
POSTAGE WILL BE PAID BY ADDRESSEE

Mercury Medical  
P.O. Box 20000  
St Petersburg FL 33742-9851



# ONLY ONE

Introducing the  
**DISPOSABLE**

## FLOWSAFE II

DISPOSABLE CPAP SYSTEM

Ideal for  
Post Anesthesia Care

Safe. Sure. Superior. CPAP

A REVOLUTION IN PACU CARE

# FLOWSAFE II

DISPOSABLE CPAP SYSTEM

- Consistent CPAP Pressure
- Disposable - No Capital Equipment Costs
- CO<sub>2</sub> Sampling Tee Fits In-Line
- No Cleaning
- Disposable Color-Coded Manometer - Verifies Delivered Pressure

# ONLY ONE

What if a CPAP device became ...

DISPOSABLE

FASTER

SAFER

MORE COMFORTABLE

MORE CONVENIENT &

MORE EFFICIENT



DISPOSABLE CPAP SYSTEM  
**FLOWSAFE II**

A REVOLUTION IN POST ANESTHESIA CARE

# ONLY ONE

11300 - 49th Street North  
Clearwater, FL 33762-4807

*Mercury Medical*

PRST STD  
U.S. POSTAGE  
PAID  
PERMIT #579  
ST. PETERSBURG, FL



# ONLY ONE

## EXPERIENCE THE ADVANTAGES!

Learn more at [www.mercurymed.com](http://www.mercurymed.com)

- View our video
- Request a free sample



**Mercury Medical**<sup>®</sup>

Yes, I would like more information about the advantages of the new Disposable Flow-Safe<sup>®</sup> II CPAP System.

Name \_\_\_\_\_

Company \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_

Zip \_\_\_\_\_ Phone \_\_\_\_\_

Email \_\_\_\_\_



**ADVANTAGE N° 1:  
BUILT-IN  
MANOMETER  
& PRESSURE RELIEF  
VALVE**  
Verifies delivered  
CPAP Pressure.



**ADVANTAGE N° 2:  
NEBULIZER (IN-LINE CAPABILITY)**  
Clinicians can administer meds without the need for mask removal.  
(Supplied Separately)



**ADVANTAGE N° 3:  
CAPNOGRAPHY (IN-LINE CAPABILITY)**  
CO<sub>2</sub> Sampling Port fits In-Line.  
(Supplied Separately)



**ADVANTAGE N° 4:  
ADVANCED MASK DESIGN**  
Lightweight contoured mask and elastic head harness provide a better seal and comfort.



**ADVANTAGE N° 5:  
USES OVER 50%  
LESS OXYGEN  
CONSUMPTION**  
With high FiO<sub>2</sub> delivery and uses standard flowmeters.

### A REVOLUTION IN POST ANESTHESIA CARE

# FLWSAFE<sup>®</sup> II

DISPOSABLE CPAP SYSTEM

NEW DESIGN

## A COMPLETE DISPOSABLE CPAP SYSTEM

